

Product Datasheet

CIITA antibody (orb387503)

Description	Rabbit polyclonal antibody to CIITA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human CIITA.
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	124
Uniprot ID	P33076
Expiration Date	12 months from date of receipt.