


Product Datasheet

COMMD4 antibody (orb387504)

Description	Rabbit polyclonal antibody to COMMD4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human COMMD4
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	22
Uniprot ID	Q9H0A8
Expiration Date	12 months from date of receipt.



Western blot of total cell extracts from...