

Product Datasheet

IL12A antibody (orb387508)

Description	Rabbit polyclonal antibody to IL12A
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of mouse IL12A
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	24
Uniprot ID	P43431
Expiration Date	12 months from date of receipt.