

Product Datasheet

SUGT1 antibody (orb387525)

Description	Rabbit polyclonal antibody to SUGT1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human SUGT1
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	41
Uniprot ID	Q9Y2Z0
Expiration Date	12 months from date of receipt.