

Product Datasheet

ARL6IP1 antibody (orb387532)

Description	Rabbit polyclonal antibody to ARL6IP1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human ARL6IP1
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	23
Uniprot ID	Q15041
Expiration Date	12 months from date of receipt.