

Product Datasheet

TTL antibody (orb387541)

Description	Rabbit polyclonal antibody to TTL
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human TTL
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	43
Uniprot ID	Q8NG68
Expiration Date	12 months from date of receipt.