

Product Datasheet

ROR2 antibody (orb387557)

Description	Rabbit polyclonal antibody to ROR2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human ROR2
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	105
Uniprot ID	Q01974
Expiration Date	12 months from date of receipt.