

Product Datasheet

RAD51B antibody (orb387561)

Description	Rabbit polyclonal antibody to RAD51B
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human RAD51B
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	42
Uniprot ID	O15315
Expiration Date	12 months from date of receipt.