

Product Datasheet

SF3B1 antibody (orb387575)

Description	Rabbit polyclonal antibody to SF3B1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human SF3B1
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	146
Uniprot ID	075533
Expiration Date	12 months from date of receipt.