

Product Datasheet

ZRSR2 antibody (orb387579)

Description	Rabbit polyclonal antibody to ZRSR2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human ZRSR2
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	58
Uniprot ID	Q15696
Expiration Date	12 months from date of receipt.