

Product Datasheet

ECSIT antibody (orb387599)

Description	Rabbit polyclonal antibody to ECSIT
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human ECSIT
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	49
Uniprot ID	Q9BQ95
Expiration Date	12 months from date of receipt.