

Product Datasheet

RNF34 antibody (orb387608)

Description	Rabbit polyclonal antibody to RNF34
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of mouse RNF34
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	42
Uniprot ID	Q99KR6
Expiration Date	12 months from date of receipt.