

Product Datasheet

HEY2 antibody (orb387611)

Description	Rabbit polyclonal antibody to HEY2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human HEY2
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	36
Uniprot ID	Q9UBP5
Expiration Date	12 months from date of receipt.