

Product Datasheet

BUB1B antibody (orb387635)

Description	Rabbit polyclonal antibody to BUB1B
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human BUB1B
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	120
Uniprot ID	O60566
Expiration Date	12 months from date of receipt.