

Product Datasheet

CD01 antibody (orb387648)

Description	Rabbit polyclonal antibody to CD01
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human CD01
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	23
Uniprot ID	Q16878
Expiration Date	12 months from date of receipt.