

Product Datasheet

PDS5A antibody (orb387657)

Description	Rabbit polyclonal antibody to PDS5A
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of mouse PDS5A
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	151
Uniprot ID	Q6A026
Expiration Date	12 months from date of receipt.