

Product Datasheet

SRSF9 antibody (orb387658)

Description	Rabbit polyclonal antibody to SRSF9
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human SRSF9
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	26
Uniprot ID	Q13242
Expiration Date	12 months from date of receipt.