

Product Datasheet

DLX5 antibody (orb387664)

Description	Rabbit polyclonal antibody to DLX5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human DLX5
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	32
Uniprot ID	P56178
Expiration Date	12 months from date of receipt.