

Product Datasheet

MED27 antibody (orb387668)

Description	Rabbit polyclonal antibody to MED27
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human MED27
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	35
Uniprot ID	Q6P2C8
Expiration Date	12 months from date of receipt.