

Product Datasheet

QPCT antibody (orb387683)

Description	Rabbit polyclonal antibody to QPCT
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human QPCT
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	38
Uniprot ID	Q16769
Expiration Date	12 months from date of receipt.