

Product Datasheet

FOXL2 antibody (orb387695)

Description	Rabbit polyclonal antibody to FOXL2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human FOXL2
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	39
Uniprot ID	P58012
Expiration Date	12 months from date of receipt.