

Product Datasheet

METTL3 antibody (orb387710)

Description	Rabbit polyclonal antibody to METTL3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human METTL3
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	64
Uniprot ID	Q86U44
Expiration Date	12 months from date of receipt.