

Product Datasheet

TPCN2 antibody (orb387711)

Description	Rabbit polyclonal antibody to TPCN2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human TPCN2
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	121
Uniprot ID	Q8NHX9
Expiration Date	12 months from date of receipt.