

Product Datasheet

FIP1L1 antibody (orb387715)

Description	Rabbit polyclonal antibody to FIP1L1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human FIP1L1
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	67
Uniprot ID	Q6UN15
Expiration Date	12 months from date of receipt.