

Product Datasheet

NOS1 antibody (orb387719)

Description	Rabbit polyclonal antibody to NOS1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human NOS1
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	161
Uniprot ID	P29475
Expiration Date	12 months from date of receipt.