

Product Datasheet

NFATC3 antibody (orb387751)

Description	Rabbit polyclonal antibody to NFATC3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human NFATC3
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	116
Uniprot ID	Q12968
Expiration Date	12 months from date of receipt.