

Product Datasheet

TNC antibody (orb387754)

Description	Rabbit polyclonal antibody to TNC
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human TNC
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	241
Uniprot ID	P24821
Expiration Date	12 months from date of receipt.