


Product Datasheet

ETV1 antibody (orb387757)

Description	Rabbit polyclonal antibody to ETV1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human ETV1
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	55
Uniprot ID	P50549
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain tissue using ETV1 antibody.