

Product Datasheet

PIRB antibody (orb387780)

Description	Rabbit polyclonal antibody to PIRB
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of mouse PIRB
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	93
Uniprot ID	P97484
Expiration Date	12 months from date of receipt.