

Product Datasheet

DAND5 antibody (orb387787)

Description	Rabbit polyclonal antibody to DAND5
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of mouse DAND5
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	20
Uniprot ID	Q76LW6
Expiration Date	12 months from date of receipt.