

Product Datasheet

PPP2R2C antibody (orb387791)

Description	Rabbit polyclonal antibody to PPP2R2C
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human PPP2R2C
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	52
Uniprot ID	Q9Y2T4
Expiration Date	12 months from date of receipt.