

Product Datasheet

LDB1 antibody (orb387793)

Description	Rabbit polyclonal antibody to LDB1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human LDB1
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	47
Uniprot ID	Q86U70
Expiration Date	12 months from date of receipt.