

Product Datasheet

ITGAV antibody (orb387808)

Description	Rabbit polyclonal antibody to ITGAV
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human ITGAV
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	116
Uniprot ID	P06756
Expiration Date	12 months from date of receipt.