

Product Datasheet

TRAK1 antibody (HRP) (orb38781)

Description	TRAK1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Trafficking kinesin-binding protein 1 protein (101-300AA)
Target	TRAK1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UPV9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.