

Product Datasheet

IFI204 antibody (orb387826)

Description	Rabbit polyclonal antibody to IFI204
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of mouse IFI204
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	72
Uniprot ID	P15092
Expiration Date	12 months from date of receipt.