

Product Datasheet

OSBPL5 antibody (orb387830)

Description	Rabbit polyclonal antibody to OSBPL5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human OSBPL5
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	99
Uniprot ID	Q9H0X9
Expiration Date	12 months from date of receipt.