

Product Datasheet

ZC3H12A antibody (orb387831)

Description	Rabbit polyclonal antibody to ZC3H12A
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human ZC3H12A
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	66
Uniprot ID	Q5D1E8
Expiration Date	12 months from date of receipt.