

Product Datasheet

ZMYND11 antibody (orb387835)

Description	Rabbit polyclonal antibody to ZMYND11
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human ZMYND11
Preservatives	0.01% NaN ₃
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	affinity purification
MW	71
Uniprot ID	Q15326
Expiration Date	12 months from date of receipt.