

## Product Datasheet

### Tubulin beta 1 antibody (orb387839)

|                            |                                              |
|----------------------------|----------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Tubulin beta 1 |
| <b>Species/Host</b>        | Rabbit                                       |
| <b>Reactivity</b>          | Human, Mouse                                 |
| <b>Conjugation</b>         | Unconjugated                                 |
| <b>Tested Applications</b> | WB                                           |
| <b>Immunogen</b>           | N-terminal region of human TUBB1             |
| <b>Preservatives</b>       | 0.01% NaN <sub>3</sub>                       |
| <b>Storage</b>             | -20°C                                        |
| <b>Note</b>                | For research use only                        |
| <b>Clonality</b>           | Polyclonal                                   |
| <b>Purity</b>              | Affinity purification                        |
| <b>MW</b>                  | 50                                           |
| <b>Uniprot ID</b>          | <b>Q9H4B7</b>                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.              |