

Product Datasheet

OPG antibody (orb387845)

Description	Rabbit polyclonal to OPG which stands for osteoprotegerin and is made up of osteoblastic cells. This...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human TNFRSF11B
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	46
Uniprot ID	O00300
Dilution Range	WB: 1:1000-3000
Expiration Date	12 months from date of receipt.