

Product Datasheet

CD133 antibody pair (orb387850)

Description	Rabbit polyclonal antibody to CD133
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal [PROM1 (N) (R2946-1)] or the C-terminal [PROM1 (C) (R2946-2)]
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	97
Uniprot ID	O43490
Expiration Date	12 months from date of receipt.