

Product Datasheet

CD31 antibody pair (orb387853)

Description	Rabbit polyclonal antibody to CD31
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of mouse CD31 or the middle region of human CD31
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	81
Uniprot ID	Q08481, P16284
Expiration Date	12 months from date of receipt.