

Product Datasheet

GPC3 antibody (orb387857)

Description	Rabbit polyclonal antibody to GPC3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human GPC3
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	65
Uniprot ID	Q9NQ94
Expiration Date	12 months from date of receipt.