

Product Datasheet

Elastin antibody pair (orb387862)

Description	Rabbit polyclonal antibody to Elastin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal or the middle region of human ELN
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	68
Uniprot ID	P15502
Expiration Date	12 months from date of receipt.