


Product Datasheet

MSH Alpha antibody (orb387875)

Description	Rabbit polyclonal antibody to MSH Alpha
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human POMC
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	29
Uniprot ID	P01189
Expiration Date	12 months from date of receipt.



Western blot of total cell extracts from...