

Product Datasheet

alpha 1 Antitrypsin antibody (orb387883)

Description	Rabbit polyclonal antibody to alpha 1 Antitrypsin
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SERPINA1
Preservatives	0.01% NaN3
Storage	-20 °C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	47
Uniprot ID	P01009
Expiration Date	12 months from date of receipt.