



Product Datasheet

CD137 antibody (orb387887)

Description	Mouse monoclonal antibody to CD137
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FC, IP
Immunogen	Recombinant human CD137 ectodomain
Target	CD137
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	Low Endotoxin
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q07011
Entrez	3604
Dilution Range	Functional application: Inhibition of cytokine production. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of cells stained using anti-h...



Flow cytometry surface staining pattern ...