


Product Datasheet

Cyclin D1 antibody (orb387896)

Description	Mouse monoclonal antibody to Cyclin D1
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, IHC-P, IP, WB
Immunogen	recombinant human cyclin D1 (amino acids 1-295)
Target	Cyclin D1
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry: Recommended dilution: 5-10 µg/ml. Flow cytometry: Recommended dilution: 3-12 µg/ml. Intracellular staining.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P24385
Entrez	595
Dilution Range	Immunohistochemistry: Recommended dilution: 5-10 µg/ml. Flow cytometry: Recommended dilution: 3-12 µg/ml. Intracellular staining.
Expiration Date	12 months from date of receipt.



Separation of HeLa cells stained using a...