


Product Datasheet


Cyclin D1 antibody (PE) (orb387897)

Description	Mouse monoclonal antibody against Cyclin D1 conjugated to PE
Reactivity	Human, Mouse, Primate, Rat
Conjugation	PE
Tested Applications	FC
Immunogen	recombinant human cyclin D1 (amino acids 1-295)
Target	Cyclin D1
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.1 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1 µg/ml. Intracellular staining.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P24385
Entrez	595
Dilution Range	Flow cytometry: Recommended dilution: 1 µg/ml. Intracellular staining.
Expiration Date	12 months from date of receipt.



Separation of HUVEC cells stained using

...



Flow cytometry intracellular staining pa...