


Product Datasheet

CD87 antibody (orb387906)

Description	Mouse monoclonal antibody to CD87
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	human myeloid cell line THP-1
Target	CD87
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q03405
Entrez	5329
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...